

UNIVERSITY
OF
CALIFORNIA

First
For UC

Now through December 31st

Special pricing on Thermo Scientific TSX ultra-low freezers

Purchase a new Thermo Scientific™ TSX600 ultra-low freezer for \$11,500 or the TSX400 for \$9,968. In addition, for each rack you purchase, get a second rack for free up to the freezer's capacity.*
A combined savings of up to \$3,452 and \$4,525!

- Low energy usage: 8.7 kWh/day for the TSX600 and 7.9 kWh/day for the TSX400
- SNAP compliant
- Unique V-drive technology is designed to reduce energy by up to 50% when compared to conventional refrigerant ultra-low freezers without compromising sample protection
- Natural refrigerants for lower environmental impact and higher cooling efficiency
- Water-blown foam insulation eliminates the refrigerant out-gassing found in other foam freezers
- Whisper quiet operation – up to 20x quieter than previous generation ULTs
- All TSX freezers are manufactured in an award-winning, zero waste to landfill facility

* Offer is valid in the United States only. August 1 through December 31, 2016. Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).



Contact your Fisher Scientific representative today!



Specifications and ordering information for Thermo Scientific TSX Series Upright Freezers (-50°C to -86°C)

Specifications

Cat. No.	Model No.	Electrical	Amps/breaker (plug)	Interior dimensions		Exterior dimensions		Shipping weight lbs. (kg)
				H x D x W in. (cm)	H x D x W in. (cm)			
NC1159161	TSX400A	115V 60Hz	8.5/15 (NEMA 5-15)	51.2 x 27 x 23.2 (130 x 68.6 x 58.9)	78.0 x 37.8 x 32.2 (198.1 x 96 x 81.8)	760 (354)		
NC1158384	TSX400D	208-230V 50/60Hz	4/10 (NEMA 6-15)	51.2 x 27 x 23.2 (130 x 68.6 x 58.9)	78.0 x 37.8 x 32.2 (198.1 x 96 x 81.8)	730 (332)		
NC1159162	TSX600A	115V 60Hz	8.5/15 (NEMA 5-15)	51.2 x 27 x 34.4 (130 x 68.6 x 87.4)	78.0 x 37.8 x 43.4 (198.1 x 96 x 110.2*)	884 (401)		
NC1158385	TSX600D	208-230V 50/60Hz	4/10 (NEMA 6-15)	51.2 x 27 x 34.4 (130 x 68.6 x 87.4)	78.0 x 37.8 x 43.4 (198.1 x 96 x 110.2*)	854 (388)		

* Door opening clearance is 34.5" (86cm)

Racks For each rack you purchase, get a second rack free up to freezer's capacity

Cat. No. Description	Dimensions H x W x D in. (cm)	Storage	TSX400 TSX600	
			TSX400	TSX600
2" box racks				
NC1158387		Boxes per rack	25	25
Adjustable Side Access Rack for 2" Boxes	11.6 x 5.4 x 26.75 (29.5 x 13.7 x 67.9)	Racks per shelf	4	6
		Racks per freezer	16	24
		Boxes per freezer	400	600

Cat. No. Description	Dimensions H x W x D in. (cm)	Storage	TSX400 TSX600	
			TSX400	TSX600
3" box racks				
NC1158386		Boxes per rack	15	15
Adjustable Side Access Rack for 3" Boxes	11.6 x 5.4 x 26.75 (29.5 x 13.7 x 67.9)	Racks per shelf	4	6
		Racks per freezer	16	24
		Boxes per freezer	240	360

Options (Field-installed must be installed by a qualified professional)

Product		Cat. No.	
		TSX400	TSX600
LN ₂ Backup System Maintains temperature down to -80°C with liquid nitrogen	Factory installed	LN4567	
	Field installed	FLN4567	
CO ₂ Backup System Maintains temperature down to -67°C with CO ₂	Factory Installed	CO4567	
	Field Installed	FCO4567	
Chart Recorder (Inkless) 6" (15.2 cm), seven day inkless recorder, -115 to +50°C, +5°C resolution	Factory Installed	CR400TSX	CR567TSX
	Field Installed	FCR400TSX	FCR567TSX
Chart Recorder (Ink) 6" (15.2 cm), seven day ink recorder, -100 to +38°C, +2°C resolution	Factory Installed	CRP400TSX	CRP567TSX
	Field Installed	CRP400TSX	FCR567TSX
Access Key Option Card access control. Includes five key cards. Supports ISO15693 and ISO14443 protocols	Factory Installed	RAC34567	
	Field Installed	FFAC34567	
Five Shelf Option	Factory Installed	5IDTSX	
	Field Installed	SK400TSX	SK600TSX
Stainless Steel Freezer Interior	Factory Installed	SS34567	

Accessories

Product	Cat. No.	
	TSX400	TSX600
Access Key Pack (North America) Includes five cards; supports ISO15693 protocol	ACU34567	
Chart Paper Ink (pack of 50)	17020	
Chart Paper Inkless (pack of 50)	6185	
Replacement Air Filter	AF34567	
Replacement Back-Up Battery	400159	
Alarm Delay Module	6903	
Cryo Gloves™ Medium	4425	
Cryo Gloves Large	4426	
Seismic Restraint Kit	TF-ULT400	TF-ULT600
Twist Lock Plugs		
Twist Lock Plug (NEMA L6-15P for 208-230V, 60Hz)	Factory Installed	US230V15ATL
	Field Installed	TL230V15A
Twist Lock Plug (NEMA L5-20P for 115V, 60Hz)	Factory Installed	US125V20ATL
	Field Installed	TL125V20A

Inquire about other Thermo Scientific freezers including Revco ExF and Forma 900 Series.

Find out more at fishersci.com



In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: fishersci.com

For Research Use Only. Not for use in diagnostic procedures. © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **0816**

